


IFE 2 x 3

Intra-Fascicular nerve Electrodes (IFE) are designed for recording peripheral nerve signals. Depending on the electrode-nerve coupling, mass signals or single spikes can be recorded.

**General Information:**

- Learn more about our materials on page 38.
- For more information about connectivity options (cables, connectors etc.), see page 31.
- Product can be supplied sterile.

Product Configuration	Catalogue Number	Connectivity Options
Platinum	9251.1009	<p>Several connection systems available, e.g.</p> <ul style="list-style-type: none"> • Omnetics – on round cable • ZIF connector – on round cable <p>Other connectivity options are available upon request. Please see page 31 for further information.</p>

✉ sales@cortec-neuro.com
 🌐 www.cortec-neuro.com
 ☎ +49 761 897 69 45 20 / 21

Copyright on all information and drawings is held by CorTec GmbH. Drawings represent the product in an abstract way and are not necessarily to scale. **All dimensions in mm.** For pricing, please request a quotation.