


Matrix Cuff 7 x 8

Polyimide technology allows the manufacturing of cuff electrodes as well. The following matrix design contains 56 electrode contacts that completely surround the nerve and that can be used for recording and stimulation.



R: reference electrode

View:
Contact openings
are shown face up

- General Information:**
- When wrapped up, anodes result in 7 rows, each with 8 individually accessible electrode contacts arranged around the circumference of the nerve.
 - Learn more about our materials on page 38.
 - For more information about connectivity options (cables, connectors etc.), see page 31.

Product Configuration	Catalogue Number
Platinum	9251.1021
Platinum Black	9251.1022

Connectivity Options
For this design, we recommend: <ul style="list-style-type: none">• D-Connector 37-way (2x) – on PCB board Please see page 31 for further information.

-  sales@cortec-neuro.com
-  www.cortec-neuro.com
-  +49 761 897 69 45 20 / 21

Copyright on all information and drawings is held by CorTec GmbH. Drawings represent the product in an abstract way and are not necessarily to scale. **All dimensions in mm.** For pricing, please request a quotation.