



ARRANGEMENT OF 360° ELECTRODE CONTACTS ON FLEXIBLE OR PUSHABLE TUBES FOR SUBCUTANEOUS OR SPINAL CORD RECORDING AND STIMULATION. PLEASE CONTACT US WITH DESIGN REQUESTS!

Example of a Tubular Electrode with 4 ring-shaped contacts.



TUBULAR ELECTRODES

PRODUCT CONFIGURATION

Please contact us with your specific design request!

CONNECTION OPTIONS

Optimal interconnection for this product will be specified according to the individual layout and application.



Copyright on all information and drawings is held by CorTec GmbH | Drawings represent abstractions of the product and are not necessarily to scale | All dimensions in mm | For pricing please request a quotation.