



Fetz Spinal Cord 8

Strip electrode for recording from and stimulation of the spinal cord, developed in cooperation with Eberhard Fetz. Please contact us for other electrode configurations or designs.

**View:**

Contact openings are shown face up

General Information:

- Openings at shaft and tip facilitate fixation.
- Product is supplied with round cable incl. strain relief.
- Learn more about our materials on page 38.
- For more information about connectivity options (cables, connectors etc.), see page 18f.
- Product can be supplied sterile.

Product Configuration**Catalogue Number**

Platinum-Iridium

1031.2024

Connectivity Options

Several connection systems available, e.g.

- Omnetics
- Touch-proof plugs
- Open cable ends

Please see page 18f. for further information.



sales@cortec-neuro.com



www.cortec-neuro.com



+49 761 897 69 45 20 / 21

Copyright on all information and drawings is held by CorTec GmbH. Drawings represent the product in an abstract way and are not necessarily to scale. **All dimensions in mm.** For pricing, please request a quotation.