


MEMS Grid Array 252

The MEMS Grid Array 252 is designed to cover a complete hemisphere with a total of 252 channels. The finger structure adjusts to the curvature of the cortex.

**View:**

Contact openings
are shown face up

General Information:

- Learn more about our materials on page 38.
- For more information about connectivity options (cables, connectors etc.), see page 31.
- Product can be supplied sterile.

Product Configuration**Catalogue Number****Platinum**

9251.1001

Product can be coated with Iridium Oxide upon request.

Connectivity Options

For this design, we recommend:

- Omnetics – direct to grid

Please see page 31 for further information.



sales@cortec-neuro.com



www.cortec-neuro.com



+49 761 897 69 45 20 / 21

Copyright on all information and drawings is held by CorTec GmbH. Drawings represent the product in an abstract way and are not necessarily to scale. **All dimensions in mm.** For pricing, please request a quotation.