


Fries-Lewis FlexShaft 12

The next generation of shaft electrodes for spike recordings: 12 channels on a flexible and thin depth electrode. As insertion aids, two types of inserters developed and designed by Ernst Strüngmann Institute are available.

**General Information:**

- Learn more about our materials on page 38.
- For more information about connectivity options (cables, connectors etc.), see page 31.
- Product can be supplied sterile.

View:

Contact openings are shown face up

The product will be supplied with a FlexShaft Inserter developed and designed by Ernst Strüngmann Institute (patent application pending).

Product Configuration	Catalogue Number
Platinum	9251.1005

Connectivity Options

For this design, we recommend:

- Harwin M50-310 connector

Other connectivity options are available upon request.
Please see page 31 for further information.



sales@cortec-neuro.com



www.cortec-neuro.com



+49 761 897 69 45 20 / 21

Copyright on all information and drawings is held by CorTec GmbH. Drawings represent the product in an abstract way and are not necessarily to scale. **All dimensions in mm.**
For pricing, please request a quotation.