


Micro Cuff 2 x 4 + 2




View:
Contact openings
are shown face up

Polyimide technology allows the manufacturing of cuff electrodes as well. The following design can be used for recording and stimulation. It will be wrapped around the nerve with a diameter of 200 μm .

**General Information:**

- When wrapped up, this cuff design contains two rows of four electrode contacts plus two ring electrodes of a diameter of 200 μm .
- Learn more about our materials on page 38.
- For more information about connectivity options (cables, connectors etc.), see page 31.

**Product Configuration****Catalogue Number**

Platinum

9251.1011

Connectivity Options

Several connection systems available, e.g.

- Omnetics - on round cable
- ZIF connector - on round cable

Please see page 31 for further information.



sales@cortec-neuro.com



www.cortec-neuro.com



+49 761 897 69 45 20 / 21

Copyright on all information and drawings is held by CorTec GmbH. Drawings represent the product in an abstract way and are not necessarily to scale. **All dimensions in mm.** For pricing, please request a quotation.