


Recording Cuff 8 + 2

Polyimide technology allows the manufacturing of cuff electrodes as well. The following designs are optimized for recording. They will be wrapped around the nerve. Varying dimensions are adapted to the peripheral nerve anatomy.


Recording Cuff 8 + 2 | Ø 1 mm



R: reference electrode
A: anodes
C: cathodes

View:
Contact openings
are shown face up

Recording Cuff 8 + 2 | Ø 2 mm

**General Information:**

- When wrapped up, anodes result in two ring contacts plus one row of single cathode contacts.
- Learn more about our materials on page 38.
- For more information about connectivity options (cables, connectors etc.), see page 31.

Product Configuration	Catalogue Number	
Recording Cuff 8 + 2 Ø 1 mm	Platinum	9251.1017
	Platinum Black	9251.1018
Recording Cuff 8 + 2 Ø 2 mm	Platinum	9251.1019
	Platinum Black	9251.1020

Connectivity Options

Several connection systems available, connected via round cable.

Please see page 31 for further information.

✉ sales@cortec-neuro.com
🌐 www.cortec-neuro.com
☎ +49 761 897 69 45 20 / 21

Copyright on all information and drawings is held by CorTec GmbH. Drawings represent the product in an abstract way and are not necessarily to scale. **All dimensions in mm.** For pricing, please request a quotation.