



Micro 8 Hexagonal
Micro 16 Hexagonal
Micro 32 Hexagonal


°AirRay research in the smallest version of ECoG grids in hexagonal arrangement. Please contact us for other electrode configurations or designs.



View:
Contact openings
are shown face up



View:
Contact openings
are shown face up



View:
Contact openings
are shown face up

General Information:

- Learn more about our materials on page 38.
- For more information about connectivity options (cables, connectors etc.), see page 18f.

Product Configuration		Catalogue Number
Micro 8 Hexagonal	Platinum-Iridium	1031.2006
Micro 16 Hexagonal	Platinum-Iridium	1031.2008
Micro 32 Hexagonal	Platinum-Iridium	1031.2028

Connectivity Options

For these products, we recommend the following connectivity option:

- Omnetics - on PCB board

Please see page 18 for further information.

Connectable to MCS Wireless Headstages.
See page 19 for more details.



sales@cortec-neuro.com



www.cortec-neuro.com



+49 761 897 69 45 20 / 21

Copyright on all information and drawings is held by CorTec GmbH. Drawings represent the product in an abstract way and are not necessarily to scale. **All dimensions in mm.** For pricing, please request a quotation.